

proFLU™ plus

Product available as:

CE Marked (*outside US*)

FDA 510(k) Cleared (*in US*)

A Multiplex Real Time PCR Assay for Detection of

- Influenza A (Matrix Gene)
- Influenza B (Non-Structural Gene)
- RSV (Polymerase Gene)

Easy to Use

- Includes Internal Control
- Includes Enzymes
- Royalty-Free
- Simple Workflow

A Superior Solution

- Faster and easier than culture
 - A few hours, not days
 - Highly automated
- Cost-Effective versus rapid tests

Clinical Performance

		Viral Culture		
		Positive	Negative	Total
ProFlu+	Influenza A Positive	132	55	187
	Influenza A Negative	0	699	699
	Total	132	754	886

Please see package insert for detailed information (Influenza B and RSV performance, discrepant analysis and other information).

Prodesse offers this product with patent rights licensed from various sources. Please see the package insert for full details.

(262) 446-0700
(888) 589-6974
Fax: (262) 446-0600
W229 N1870 Westwood Drive
Waukesha, WI 53186
www.prodesse.com



Prodesse

Real Time Solutions